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**NMIJ CRM is distributed from  
National Metrology Institute of Japan (NMIJ).**

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#### **Features of NMIJ CRM**

- Production in accordance with quality management system based on ISO 17034
- Certification in accordance with ISO Guide 35
- Determination of analyte by state of the art technology including primary methods
- Traceable to SI
- Uncertainty evaluation in accordance with Guide to the Expression of Uncertainty in Measurement (GUM)

#### **NMIJ CRM fits for**

- Calibration of instruments
- Validation of analytical methods and measurement results
- Quality Control for analytical methods.

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Reference material characterized by a metrologically valid procedure for one or more specified properties, accompanied by a reference material certificate that provides the value of the specified property, its associated uncertainty, and a statement of the metrological traceability

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Material, sufficiently homogeneous and stable with respect to one or more specified properties, which has been established to be fit for its intended use in a measurement process

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## Utilization of Reference Materials

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Reference materials (RMs) are generally used for ;

- 1) Calibration of analytical/measuring equipment,
- 2) Assignment of chemical/physical quantities of materials,
- 3) Evaluation of analytical/measuring methods, and
- 4) Confirmation of skill of analysts/performance of organizations.

1) “Calibration” means adjustment of a value indicated by the analytical/measuring equipment to a correct (or true) value. Calibration also includes conversion of an indicated value with a physical quantity into concentration or another quantities by using a calibration curve.

2) “Assignment” is an analysis of a sample in comparison with a certified value of a certified reference material (CRM) or with that of an indicative value of an RM, and so includes the second meaning of “Calibration”. For example, in the case of titrimetric analysis, a concentration of unknown sample is determined on the basis of an RM whose purity has been certified. The measurement law of Japan defines a “Calibration” as assignment of “lower (in the metrological traceability system)” RMs by using “higher” RMs. On the other hand, high-purity and other similar CRMs are usually used for “Calibration” and “Assignment”.

3) “Evaluation” is a judgment of whether the analyst’s analytical/measuring method is appropriate and reliable or not. “Evaluation” may often be called “Validation”.

In the case of 4) “Confirmation”, it is confirmed whether an organization/analyst has sufficient technical skills to obtain reliable experimental data. In general, matrix CRMs including environmental, nutritional, ..., etc., are used for the purposes of “Evaluation” and “Confirmation”.

A CRM is an RM with an RM certificate, that is described in ISO Guide 30: Reference materials - Selected terms and definitions. The certificate includes certified value(s) and their uncertainty(s), that have been determined by procedures to ensure traceability to the International System of Units (SI). Thus, a CRM is absolutely necessary for an evaluation of the accuracy of an analytical value and for confirmation of traceability to SI units. In addition, in an accreditation system based on ISO/IEC 17025, CRMs are useful for the confirmation of traceability to SI units.

If CRMs are used inappropriately, their merits will not be recognized. In general, CRMs have good homogeneity and stability. However, the long-term stability differs between individual CRMs, and some CRMs have a short period of validity. Please pay attention to the expiry date and storage conditions described in the RM certificate provided with each CRM. For several CRMs, drying and storage rules before using are strictly described in the RM certificates.

The RM certificate of each CRM includes the scope, certified value, sample form, analytical method used for certification, preparation method, instructions for use, and so on. In some instances, indicative values are included in the certificates. The indicative value is the property value which does not meet the NMIJ’s requirements of the certified value including a case that it was estimated that evaluations of the analytical method used for determination of property value are not enough. NMIJ provides the RMs as NMIJ RM which is assigned only indicative value.

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The mole (unit: mol) was adopted as an SI base unit by the 14th General Conference on Weights and Measures (1971) (14<sup>th</sup> Conference Generale des Poids et Mesures), convened under the Meter Convention. Traceability to amount of substance can be achieved through continuous (unbroken) chains of SI to the CRM (the primary reference material), and those of the CRM to other RMs (by using analytical equipment, etc.). The development of the primary reference material, which is directly linked to an SI base unit, requires diligence and a high level of skill. In many cases, the primary reference materials have been developed by National Metrology Institute (NMI).

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## Method for the Determination of Property Values of CRMs

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One method to determine property values is a primary method of measurement. A primary method is defined as “a method having the highest metrological properties, whose operation can be completely described and understood, for which a complete uncertainty statement can be written down in terms of SI units”. The primary methods identified by the Consultative Committee for Amount of Substance (CCQM) are coulometry, gravimetry (including gravimetric preparation method), titrimetry, isotope dilution mass spectrometry, and the freezing point (depression) method. Most of NMIJ CRMs can be accepted internationally as representative RMs traceable to SI because the property values of them are determined using these primary methods.

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# International Comparison and Global Mutual Recognition Arrangement

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Several working groups under the CCQM have been established. Each working group (WG) plans and performs various international comparisons. The Calibration and Measurement Capability (CMC) of each National Metrology Institute (NMI) and the comparability between primary standards of NMIs have been checked through the international comparisons. In the case of “reference material”, the property value is generally determined using analytical results. Therefore, in international comparisons, a pilot laboratory distributes identical samples to each NMI, and then the CMC to determine the property value of the sample is evaluated for each NMI. In such international comparison, each NMI selects analytical techniques, which means that the CRMs produced by the techniques are also evaluated simultaneously. NMIJ participates in this international comparison. In addition, NMIJ is reviewed the technical capability by the expert of the overseas National Metrology Institute (Peer review) regularly. The CMCs and the range of determining the certified value for CRM of NMIJ which are approved based on these results are registered in the key comparison database (KCDB: CIPM MRA, <https://www.bipm.org/kcdb/>) which International Bureau of Weights and Measures (BIPM) manages, and they are published on the website of BIPM. The most of NMIJ CRMs are registered in this database and internationally recognized.

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## List of Distributors and Classification

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認証標準物質 (NMIJ CRM) は、下記の事業者が取り扱っています。  
NMIJ CRMs are distributed only through the following companies.  
Please contact them to purchase.

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### Customer Registration of NMIJ CRM & RM

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<https://unit.aist.go.jp/qualmanmet/refmate/>

Completing customer registration on the following NMIJ Website will facilitate notification of any revision of the information about registered CRM(s) or RM(s).

Please register an individual CRM or RM you purchased.

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#### 【Disclosure of Personal Information to Third Parties】

NMIJ shall handle the customers personal information strictly in compliance with “Act on Protection of Personal Information” and use them for the purposes only to send notification on NMIJ CRMs and to research CRM market.

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A value after the sign “±” in each certified value is the expanded uncertainty interval calculated using a coverage factor ( $k$ ) of 2, which gives a level of confidence of approximately 95 %.

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Please check an expiration of each certification before purchase. The latest information, for example, the distribution of new CRMs, the end of distribution and lot updates, etc., and sample of certificate for each CRM are available from NMIJ Website (<https://unit.aist.go.jp/nmij/english/refmate/>).

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Electrolytic Conductivity Standard Solution  
- Aqueous Solution of Potassium Chloride  
(0.01 mol kg<sup>-1</sup>)

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CRMs 1001-a ~ 1005-a, 1006-a ~ 1010-a, 1017-a and 1018-a should be kept in dry and clean atmosphere such as a desiccator at a temperature between 5 °C and 35 °C.

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After rinsing the alloy surface with acetone or the like, CRM 1019-a and CRM 1020-a should be kept in dry and clean atmosphere such as a desiccator at room temperature.

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CRM 5203-a should be kept in a dry and clean environment at a temperature between 5 °C and 35 °C. The storage under the nitrogen flow is recommended.

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CRM 5205-a should be stored in a clean and dry environment, such as in a desiccator, at a temperature between 5 °C and 35 °C.

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
CRM 5206-a should be stored in a clean and dry environment, such as in a desiccator, at a temperature between 5 °C and 35 °C.

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CRM 5207-a should be stored at a temperature between 5 °C and 35 °C in a nitrogen atmosphere in which the relative humidity is kept below 10 %.

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RM 8158-a should be kept in its aluminum laminated plastic bag at a temperature between 15 °C and 25 °C and shielded from light.

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CRM 5601-a must be stored in a clean environment at room temperature. It is recommended to store this CRM under an atmosphere of dry air or nitrogen gas away from any radioactive sources.

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CRMs 5602-a and 5606-a should be kept light-shielded in an atmosphere of dry air or nitrogen gas at a temperature between 15 °C and 35 °C and keep it away from any radioactive sources.

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RM 5607-a should be stored at a temperature between 15 °C and 35 °C in clean environment and away from any radioactive sources.

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The CRMs above should be stored under clean environment at a temperature between 4 °C and 30 °C. Do not freeze them.

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Example of specific volume of absorbed nitrogen for each relative pressure (CRM 5714-a)

物質名 : カーボンブラック (窒素吸着量 -BET100)  
 Description : Carbon Black (Nitrogen Specific Volume Adsorbed - BET100)

相対圧 Relative Pressure (Pa/Pa)	窒素吸着量 Specific Volume of Absorbed Nitrogen (cm <sup>3</sup> /g)	相対圧 Relative Pressure (Pa/Pa)	窒素吸着量 Specific Volume of Absorbed Nitrogen (cm <sup>3</sup> /g)	相対圧 Relative Pressure (Pa/Pa)	窒素吸着量 Specific Volume of Absorbed Nitrogen (cm <sup>3</sup> /g)
0.050	24.4 ± 1.3	0.250	32.9 ± 1.0	0.600	50.2 ± 2.5
0.075	25.63 ± 0.53	0.275	34.0 ± 1.1	0.650	53.6 ± 2.8
0.100	26.70 ± 0.59	0.300	35.2 ± 1.2	0.700	57.6 ± 3.1
0.125	27.73 ± 0.66	0.350	37.4 ± 1.4	0.750	62.8 ± 3.5
0.150	28.74 ± 0.72	0.400	39.7 ± 1.6	0.800	70.0 ± 4.2
0.175	29.76 ± 0.79	0.450	42.1 ± 1.8	0.850	81.0 ± 5.5
0.200	30.79 ± 0.87	0.500	44.6 ± 2.0	0.900	105 ± 11
0.225	31.85 ± 0.94	0.550	47.2 ± 2.2		

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The RMs above show a weak photocatalytic activity. The RMs should be stored in a clean environment at a temperature between 5 °C and 35 °C with the lid of tightly closed and shielded from direct sunlight.

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Technical informations are necessary when RMs above are used as a representative test material (ISO/TS 16195:2013) in developing and validating test methods using nano-objects (for example, toxicology tests).

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The CRMs 3009-a, 3011-a and 3013-a should be stored at a temperature between 15 °C and 30 °C, at a relative humidity of 60 % or less and shielded from sunlight.

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The CRM 3012-a should be stored at a temperature between 15 °C and 35 °C and at a relative humidity of 60 % or less in a clean place shielded from light. It should not be affected by acids, bases, oxidants, reductants, organic substances and others.

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CRMs above should be kept in the glass bottle sealed in plastic bag, and should be stored in a clean environment at a temperature between 15 °C and 30 °C.

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Electrolytic Conductivity  
Standard Solution  
– Aqueous Solution  
of Potassium Chloride  
(0.01 mol kg<sup>-1</sup>)

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The solution of above CRM should be kept in the high-density polyethylene bottle sealed in an aluminum-laminated plastic bag. The CRM should be stored in a clean environment at a temperature of 10 °C or less where the solution could not be frozen.

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CRMs 4003-b and 4039-a should be stored in a clean and dark environment at a temperature between 2 °C and 10 °C.

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Compliance with Japanese "Poisonous and Deleterious Substances Control Law" is necessary because the CRMs are "Deleterious Substances" in the law.

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CRM 4001-b should be stored under clean and dark environment at a temperature of approximately -20 °C.

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CRM 4036-b and 4074-a should be stored under clean and dark environment at a temperature between -15 °C and -25 °C.

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CRMs 4055-a and 4057-a should be stored under clean and dark environment at a temperature between 2 °C and 8 °C.

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CRM 4056-a should be stored under clean and dark environment at a temperature between 0 °C and 6 °C .

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CRM 4058-a should be stored under clean and dark environment at a temperature between 2 °C and 10 °C.

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<sup>14</sup>C contents are valued for CRM 4001-b as indicative values.

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CRM 4055-a includes *tert*-butylcatechol as stabilizing material.

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Your signature to a pledge card is necessary to purchase the CRM or RM.

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Compliance with Japanese "Poisonous and Deleterious Substances Control Law" is necessary because the CRM is "Deleterious Substances" in the law.

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CRM 5401-a should be stored in a clean and dark environment at a temperature between 5 °C and 35 °C.

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CRM 4601-c should be stored at a temperature between 15 °C and 25 °C and shielded from light in a clean desiccator.

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CRM 4602-a should be stored in a clean and dark environment at a temperature between 2 °C and 10 °C.

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CRM 4603-a should be stored at a temperature between 15 °C and 35 °C and at a relative humidity of 60 % or less in a clean place shielded from light.

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CRMs 4213-a and 4220-a should be stored in a clean and dark environment at a temperature between 15 °C and 25 °C.

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CRM 4222-e should be stored in a clean and dark environment at a temperature between 2 °C and 10 °C.

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CRMs 4228-b and 4229-a should be stored in a clean and dark environment at a temperature between 15 °C and 30 °C.

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Your signature to a pledge card is necessary to purchase the CRMs.

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Compliance with Japanese “Poisonous and Deleterious Substances Control Law” is necessary because the CRM is “Deleterious Substances “ in the law.

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CRM 4213-a will end its distribution on March 31, 2025.

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Density is indicated as information.

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Potassium Hydrogen Phthalate  
for Quantitative NMR (<sup>1</sup>H)

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CRM 6001-a should be stored under dark environment at a temperature between 0 °C and 6 °C .

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CRMs 6002-a, 6003-a, and 6004-a should be stored under dark environment at a temperature between 2 °C and 8 °C .

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CRMs 6005-a, 6006-a and 6008-a should be stored at a temperature between 15 °C and 25 °C and shielded from light in a clean desiccator.

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CRM 6007-a should be stored in a clean and dark environment at a temperature between 15 °C and 25 °C. It is recommended to store CRM 6007-a at a temperature around 5 °C, if storage for a longer period is needed.

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CRMs above should be stored at a temperature between 15 °C and 25 °C and shielded from light in a clean desiccator.

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◆ NIST : National Institute of Standards and Technology (米国)

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CRM 6202-a should be stored under clean environment at a temperature of approximately 4 °C, but it should not be frozen.

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The CRMs 6204-b and 6205-b should be stored at a temperature between -30 °C and -20 °C in a freezer.

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CRM 6210-a should be stored at temperatures of ca. -80 °C and shielded from light.

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CRMs 6211-a and 6212-a should be stored at a temperature between -30 °C and -15 °C in a clean freezer and protected from light.

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
The measurement results of this RM (AIST-MAB) including the primary structure are summarized as “Technical note (case study) A” in NMIJ home page.

([https://unit.aist.go.jp/nmij/english/refmate/information/6208a\\_en\\_info.pdf](https://unit.aist.go.jp/nmij/english/refmate/information/6208a_en_info.pdf))

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CRM 5010-a should be stored in a clean and dark environment at a temperature of 5 °C or less.

\* CRM 5010-a は各重合度成分ごとの分子量、質量分率およびモル（数）分率について、認証書に記載があります。  
For CRMs 5010-a, molecular weight, mass fraction, and mole (number) fraction for each degree of polymerization are also indicated in the certificate.

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The CRMs 5011-a should be stored under clean and dark environment at a temperature of 25 °C or less. A temperature of approximately 5 °C is recommended for long-term storage.

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CRM 5010-a will end its distribution on March 31, 2025.

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RM 5012-a should be stored in a clean and dark environment at a temperature between 5 °C and 35°C.

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CRM 5012-a will end its distribution on March 31, 2025.

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CRM 7202-c should be stored in a clean and dark environment at a temperature of approximately 5 °C.

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CRMs 7204-a and 7303-a should be stored in a clean and dark environment at a temperature between 5 °C and 35 °C.

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CRM 8302-a should be stored in a clean environment at a temperature between 15 °C and 30 °C and shielded from light.

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Compliance with Japanese "Poisonous and Deleterious Substances Control Law" is necessary because the CRM is "Poisonous Substances" in the law.

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CRMs 7407-a and 7408-a should be stored in a clean environment at a temperature between -20 °C and -30 °C and shielded from light.

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CRM 7906-a should be stored in a clean and dark environment at a temperature between 15 °C and 30 °C.

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In disposing CRM7906-a, it corresponds to the low-concentration PCB waste with the Japanese legislation, Act on Special Measures concerning Promotion of Proper Treatment of PCB Wastes.

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CRM 7906-a will end its distribution on March 31, 2025.

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CRM 7901-a should be stored in a clean and dark environment at a temperature between 5 °C and 35 °C.

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CRMs 7402-a and 7405-b should be stored in a clean and dark environment at a temperature between 5 °C and 35 °C.

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CRMs 7503-b and 7533-a should be stored in a clean and dark environment at a temperature between 5 °C and 35 °C.

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CRM 7203-a should be stored in a clean environment at a temperature between 5 °C and 25 °C and shielded from light.

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CRMs 6206-b and 6207-b should be stored at temperature of -20 °C or less and shielded from light.

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CRMs 7520-a and 7521-a should be kept at a temperature between -20 °C and -30 °C and shielded from light.

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CRM 7504-a should be stored under clean and dark environment at temperature of approximately -30 °C.

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CRMs 8112-a and 8115-a should be stored in a clean and dark environment at a temperature between 15 °C and 35 °C.

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CRMs 8108-b and 8110-b should be stored in a clean and dark environment at a temperature between 5 °C and 35 °C.

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The CRM 8155-a should be stored in a clean and dark environment at a temperature between 15 °C and 25 °C.

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CRMs 8152-b and 8156-a should be stored in a clean and dark environment at a temperature between 15 °C and 35 °C.

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The CRMs above should be used in compliance with the high pressure gas safety law. They should be stored in a clean environment away from direct sunlight at a temperature between 10 °C and 35 °C. The CRMs above should be stored at a place with good ventilation and do not use fire near the cylinder. CRMs above should be secured with a chain to prevent them falling.

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Each certified value shown in table above is an example. The certified value is slightly different in each cylinder.

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All high-purity gas CRMs above are made-to-order products.

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Polyvinyl Chloride  
(Phthalate Esters in PVC Resin Pellet)  
(P.42)

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For more information, please visit our web site. <https://unit.aist.go.jp/nmi/english/refmate/>

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The CRMs above should be used in compliance with the high pressure gas safety law. They should be stored in a clean environment away from direct sunlight at a temperature between 10 °C and 35 °C. CRMs above should be secured with a chain to prevent them falling.

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All gas mixture CRMs above are made-to-order products.

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(注) 形状 1 と形状 2 は別売りです。 Shape 1 and Shape 2 are available separately.

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The RMs above are recommended to be stored in noble gas at a temperature of 23 °C ± 5 °C as well as relative humidity of 50 % or less.

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The several certified values are indicated for multiple temperatures within the operating temperature range. (see P. 47)

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Form 1 and Form 2 are available separately.

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Example of the certified values for multiple temperatures within the operating range. (CRM 5805-a)

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Description : High-purity copper for Thermal Expansivity Measurements

認証値 Certified Values		認証値 Certified Values	
温度 Temperature (K)	熱膨張率 Coefficient of Thermal Expansion (10 <sup>-6</sup> K <sup>-1</sup> )	温度 Temperature (K)	熱膨張率 Coefficient of Thermal Expansion (10 <sup>-6</sup> K <sup>-1</sup> )
20	0.279 ± 0.052	180	14.665 ± 0.029
40	2.260 ± 0.046	200	15.174 ± 0.030
60	5.457 ± 0.041	220	15.585 ± 0.030
80	8.328 ± 0.037	240	15.927 ± 0.031
100	10.486 ± 0.034	260	16.220 ± 0.032
120	12.044 ± 0.032	280	16.479 ± 0.033
140	13.182 ± 0.031	300	16.713 ± 0.034
160	14.024 ± 0.030		

熱物性標準物質 (熱膨張率)  
CRMs for Thermophysical properties  
(Coefficient of Thermal Expansion)

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CRMs 5807-a and 5809-a should be stored at a temperature of 23 °C ± 10 °C as well as relative humidity of 50 % or less.

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CRM 5810-a should be stored in a clean environment at a temperature between 5 °C and 35 °C

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The several certified values are indicated for multiple temperatures within the operating range. (see P. 47)

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CRM 5808-b should be stored in a clean environment at a temperature between 5 °C and 35 °C

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RM 1401-a is recommended to be stored at a temperature of 23 °C ± 5 °C as well as relative humidity below 50 %.

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The several Indicative values are indicated for multiple temperatures within the operating temperature range. (see P. 47)

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NMIJ RM is a reference material, but not a certified reference material.

熱物性標準物質 (熱拡散率・熱伝導率)  
CRMs for Thermophysical properties  
(Thermal Diffusivity, Thermal Conductivity)

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Electrolytic Conductivity Standard Solution  
– Aqueous Solution of Potassium Chloride (1 mol kg<sup>-1</sup>)

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Date of Shipment: XXXXX xx, 20xx

5121 b 02-240315-240315

National Institute of Advanced Industrial Science and Technology

## National Metrology Institute of Japan

## Reference Material Certificate

NMIJ CRM 5121-b02

No. +++

Electrolytic Conductivity Standard Solution  
 –Aqueous Solution of Potassium Chloride (1 mol kg<sup>-1</sup>)

This certified reference material (CRM) is produced in accordance with the NMIJ's management system and is in compliance with ISO 17034 and ISO/IEC 17025. This CRM is intended for use in the calibration of electrolytic conductivity.

**Certified Value**

The certified value of this CRM is given in the table below. The uncertainty of the certified value is the expanded uncertainty obtained by multiplying the combined standard uncertainty by a coverage factor ( $k$ ) of 2, and it is the half-width of an interval estimated to have a level of confidence of approximately 95 %.

Property	Certified value S m <sup>-1</sup>	Expanded uncertainty S m <sup>-1</sup>
Electrolytic Conductivity (25 °C)	10.867	0.060

**Analysis**

The certified value of this CRM was determined from measurements of the geometry (length and cross-sectional area) of the glass cell for electrolytic conductivity and the impedance of the solution at 25.000 °C ± 0.014 °C ( $k = 2$ ).

**Metrological Traceability**

The certified value of this CRM was determined on absolute basis; the geometry (length and cross-sectional area) of the glass cell for electrolytic conductivity was calibrated with a coordinate measuring machine calibrated as traceable value to NMIJ's national standard of length, and a LCR meter for impedance measurements was calibrated on Japan Calibration Service System (JCSS). Therefore, the certified value is traceable to the SI.

**Mutual Recognition Arrangement under Metre Convention**

The certified value of this CRM is recognized for international equivalence based on the Mutual Recognition Arrangement under the Metre Convention (CIPM MRA). The calibration measurement capability (CMC) of NMIJ related to this CRM is registered in the Key Comparison Database (KCDB) (see <https://www.bipm.org/kcdb/>) of the International Bureau of Weights and Measures (BIPM).

**Expiration of Certification**

This certificate is valid for 6 months from the date of shipment or until January 30, 2028, whichever comes earlier, provided

· The certificates of NMIJ CRMs include the information on following properties.

- Application of the CRM
- Analysis
- Indicative Value(s) (If applicable)
- Mutual Recognition Arrangement (CIPM MRA) (If applicable)
- Expiration of Certification
- Homogeneity
- Instructions for Use
- Preparation
- Certified Value(s)
- Metrological Traceability
- Description of the Material
- Instructions for Storage
- Precautions for Handling (If applicable)
- Revision history (If applicable)

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A		アルファベット		Alphabet		B	
absorbed nitrogen	17	ammonium ion	19	ammonium ion	19	B	34, 38, 39
ABS Resin	41	amyloid $\beta$	30	amyloid $\beta$	30	Ba	34, 39
acetamiprid	35	aqueous particle dispersion	16	aqueous particle dispersion	16	base	19
acrylonitrile	23	aqueous solution	20	aqueous solution	20	benzo[a]pyrene	26
Ag	34	arginine	28	arginine	28	benzyl butyl phthalate	42
AIST-MAB	32	argon	45	argon	45	BET	17
Al	13, 34, 38, 39	arsenate	38	arsenate	38	biodiesel	34
Al <sub>2</sub> O <sub>3</sub>	46	arsenic acid	37	arsenic acid	37	bis(2-ethylhexyl) adipate	42
Al <sub>2</sub> O <sub>3</sub> -TiC	48	arsenic compound	37, 38	arsenic compound	37, 38	bis(2-ethylhexyl) phthalate	42
alanine	28	arsenic-dope	13	arsenic-dope	13		42
AlAs	13	arsenite	38	arsenite	38	bisphenol A	22
albumin	29	arsenobetaine	37	arsenobetaine	37	3,5-bis(trifluoromethyl)benzoic Acid	25
aldosterone	31	arsenosugar	37	arsenosugar	37	1,4-bis(trimethylsilyl) - 2,3,5,6 - tetrafluorobenzene	25
alloy	12	artificial urine	35	artificial urine	35		
alumina ceramic	46	As	13, 34, 37, 38, 39	As	13, 34, 37, 38, 39		
amino acid sequence	32	As-dope	13	As-dope	13		
ammonium chloride	19	aspartic acid	28	aspartic acid	28		
		averaged diameter	16	averaged diameter	16		

biuret method	29	decabrominated diphenyl ether	40
BN delta-Layer	13		41
bragg angle	18	degree of polymerization	
brown rice flour	38, 39		33
butane	44	delta-layer	13
butylmethylether	23	deoxyribonucleic acid	29
		diarrhetic shellfish toxin	
<b>C</b>			40
Ca	34, 37, 38, 39	dibromochloromethane	11, 23
cadmium	39	dibutyl phthalate	42
calcium	19	1,4-dichlorobenzene	22
calcium carbonate	19	1,1-dichloroethylene	22
carbon	47	dichloromethane	22
carbon black	17	1,2-dichloropropane	22
carbon dioxide	43	dicyclohexyl phthalate	42
carbon monoxide	43	diethyl phthalate	42
carbon tetrachloride	22	diheptyl phthalate	42
CB	35	dihexyl phthalate	42
Cd	34, 37, 38, 39, 40, 41	3 $\beta$ , 4 $\alpha$ -dihydroxy-5 $\alpha$ -androstan- 17-one	30
ceramic	48	diisobutyl phthalate	42
chloride ion	19	dimethylarsinic acid	38
chlorinated paraffin	24	dimethyl phthalate	42
chlorine	24	Dinitrogen Oxide	45
chloroalkane	24	dinophysistoxin-1	36, 40
cholesterol	27	dioctyl phthalate	42
Co	34, 39	1,4-dioxane	23
cod fish	37	dipentyl phthalate	42
combustion supporting gas		dipropyl phthalate	42
	43	dissolved Silica	36
copper	46, 47	DNA	29
cortisol	31	dot-array	14
C-peptide	31	DTX1	36, 40
Cr	12, 34, 37, 38, 41	<b>E</b>	
creatinine	27	electrolytic conductivity	20
crystallite size	18	enthalpy	25, 50
Cu	34, 37, 38, 39	17 $\beta$ -estradiol	27
cyclohexane	25, 50	ethane	44
cystine	28	ethanol	23
		ethylbenzene	22
<b>D</b>			
deamidated C-peptide	31		
		etofenprox	40
		<b>F</b>	
		F	25
		fatty acid	18
		Fe	12, 34, 37, 38, 39
		Fe-Cr alloy	12
		Fe-Ni alloy	12
		fenitrothion	40
		flammable gas	43, 44
		fthalide	40
		fuel	34
		fused silica	15
		fusion enthalpy	25, 50
		<b>G</b>	
		Ga	13
		GaAs	13
		GaAs/AlAs	13
		glass-like carbon	47
		glutamic acid	28
		glycine	28
		graphite	49
		<b>H</b>	
		H	25, 26
		heavy metal	38, 41
		Hexane	45
		Hg	34, 37, 41
		Hijiki	37
		histidine	28
		homologue	24
		human serum	
			11, 29, 31, 32, 35
		hydrocortisone	27
		4-hydroxy-clomifene	30
		<b>I</b>	
		International System of Units	
			5
		isobutane	44

isobutyl group	18	molybdenum	49	Pb	21, 34, 38, 39, 41
ISO Guide 30	4, 5	monoclonal antibody	32	PCB	35
ISO/IEC 17025	4, 5	monodisperse	16	pentane	44
isoleucine	28	MTBE	23	peptide	31
isopentane	44	<b>N</b>		perfluoroalkyl substance	41
isoprothiolane	40	N <sub>2</sub> O	45	perfluorooctanesulfonic acid	41
Isotopic	21	Na	34, 37, 38, 39	perfluorooctanoic acid	23
isotopic ratio	21	neonicotinoid	35	pesticide	35, 40
isotropic graphite	49	Ni	12, 34, 37, 38, 39	phase transition enthalpy	25, 50
ISO/TS 16195	18	nickel alloy	12	phase transition temperature	25, 50
<b>K</b>		Ni-Fe alloy	12	phenylalanine	28
K	34, 37, 38, 39	NIST	29	phosphate ion	36
<b>L</b>		nitrate ion	36	phthalate ester	42, 44
lake sediment	34	nitrite ion	36	polybrominated diphenyl ether	41
latex	16	nitrogen	17, 27, 43, 45	polycarbonate	15
lead	21	nitrogen specific volume	17	polychlorinated biphenyl	35
leucine	28	adsorbed	17	poly(ethylene glycol)	33
Li	39	NMR	25, 26	poly(ethylene glycol) nonylphenyl ether	33
light scattering	16, 33	nuclear magnetic resonance	25	polystyrene	14, 16, 33, 41
liquefied gas	43, 44	number average diameter	16	polyvinyl chloride	42, 44
LNG	43	number-average molecular mass	33	positron annihilation lifetime	15
lysine monohydrochloride	28	nutrient	36	positron defect measurement	15
<b>M</b>		<b>O</b>		positron hole-size	15
MAB	32	okadaic acid	40	positron lifetime	15
mass-average molar mass	33	organic contaminant	35	potassium chloride	11, 20
mass-average molecular mass	33	ortho-positronium lifetime	15	potassium hydrogen phthalate	25, 26
melting point	25, 50	oxygen	43, 45	potassium	
methane	43, 45	<b>P</b>		perfluorooctanesulfonate	26
methanol	34	P	34, 39	progesterone	27
methionine	28	palm oil	34	proline	28
methylmercury	37	particle mass	16	protein	29
Mg	34, 37, 38, 39	particle size	16		
milk	39	particle size distribution	16		
Mn	34, 37, 38, 39				
Mo	34, 38, 39				

PVC Resin	42, 44	SRM 927d	29	<b>W</b>	
pyroglutamylated C-peptide		stainless steel	12, 15	water	26, 34, 38
	31	static light scattering	33	water suspension	16
<b>Q</b>		steroid hormone	11, 31, 32	white rice flour	38
quartz glass	48	styrene	23	<b>X</b>	
<b>R</b>		sulfur dioxide	43	xylene	22
radius of gyration	33	super lattice	13	<b>Z</b>	
Raman	14	surface area	18	zinc	19
Rb	34, 38, 39	surface modification	18	Zn	34, 37, 38, 39
ribonucleic acid	29	surface modified	18		
river water	34	<b>T</b>			
RNA	29	tap water	38		
<b>S</b>		tea leaf	39		
S	34	testosterone	27, 31		
Sb	34, 38	thiacloprid	35		
scallop edible part	40	thiamethoxam	35		
scallop midgut gland	40	threonine	28		
Se	34, 37, 38	tissue	37		
seawater	34, 36	titanium (IV) oxide	18		
seaweed	37	titanium nitride	48		
serine	28	toluene	22		
short-chain chlorinated paraffin		toxic gas	43		
	24	toxicology test	18		
Si	13	trace element	34, 37, 38, 39		
Si	4, 5, 6	trichloroacetic acid	23		
silica	36	tris (hydroxymethyl)			
silica glass	46	aminomethane	19		
silicon	15, 46, 47, 50	tungsten	14		
single crystal of silicon	46, 50	tyrosine	28	<b>U</b>	
single-crystal silicon	15, 47			unpolished rice	40
Sn	34			urea	27
solid-solid phase transition				uric acid	27
	25, 50			urine	35
soybean	39	<b>V</b>			
specific surface area	17, 18	Valine	28		
specific volume	17				
squid	39				
Sr	34, 38, 39				

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If you have any inquiries, please contact the following address.

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<https://unit.aist.go.jp/nmij>

Reference Materials Office

Center for Quality Management of Metrology

National Metrology Institute of Japan (NMIJ)

National Institute of Advanced Industrial Science and Technology (AIST)

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