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AIST Nanocharacterization Facility

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Advanced Research Infrastructure for Materials and Nanotechnology in Japan (ARIM)
by the Ministry of Education, Culture, Sports, Science and Technology (MEXT).

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AIST Nanocharacterization Facility (ANCF), one of AIST's Open Research Facilities, participates in both TIA, an open innovation hub, and the Advanced Research Infrastructure for Materials and Nanotechnology (ARIM) in Japan sponsored by the Ministry of Education, Culture, Sports, Science and Technology (MEXT). Through these two channels, we respond to a wide range of measurement requests from all over the country.

文部科学省マテリアル先端リサーチインフラ (ARIM) Advanced Research Infrastructure for Materials and Nanotechnology (ARIM) funded by MEXT

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Advanced Research Infrastructure for Materials and Nanotechnology (ARIM) is one of the research projects by Japan's Ministry of Education, Culture, Sports, Science and Technology (MEXT), which started in FY 2021 under the 10-year-long plan. ARIM is devoted to construct and maintain state-of-the-art open facilities for materials processing, measurement and analysis as well as their data collection and utilization system for public share-use, thus playing a key role in "Materials DX Platform" so as to promote Japan's materials innovations. This project utilizes nationwide network of both R&D facilities and highly skilled technical professionals, which was established under the Nanotechnology Platform Project. AIST Nanocharacterization Facility (ANCF) is a member of ARIM, together with AIST Nanoprocessing Facility (NPF), for data acquisition, analysis, and collection of cutting-edge nanometer-scale materials including quantum and electronic materials.

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TIA is an open innovation hub operated in cooperation with AIST, NIMS, the University of Tsukuba, KEK, the University of Tokyo and Tohoku University. We will comprehensively bring together our R&D capabilities (human resources, facilities, intellectual property, etc.) to provide consistent support, from knowledge creation to industrialization.

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AIST provides access to a wide range of advanced equipment and know-how for shared use, with the aim of widely disseminating such know-how and findings that have been accumulated so far.

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Development and disclosure of advanced measurement and analysis technology by AIST Nanocharacterization Facility (ANCF)

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AIST is developing advanced measurement and analysis technology with the aim of strengthening domestic industries, leading the creation of new industries, and contributing to social innovation. Our mission is to create measurement technologies that make the invisible visible, and then to apply such technologies to various fields to refine them as formed analytical methods. To achieve this mission, we are trying to solve problems in society by sharing the equipment and technology that we have developed.


At ANCF, we mainly offer our proprietary measuring equipment and technologies for common use that are not commercially available. For example, one of AIST's unique devices, the positron probe microanalyzer, uses a positron microbeam generated from an electron accelerator; there is only one other similar device in the world, in Germany. In addition, AIST has developed a superconducting detector for X-ray absorption fine structure analyzers and scanning electron microscopes to achieve higher resolution and sensitivity than existing devices. Even with transient absorption spectrometers, real surface probe microscopes, and solid-state NMR devices, we have achieved unsurpassed performance by precisely combining basic commercially-available devices and applying highly skilled sample preparation methods and analytical methods.

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Making the invisible visible
Measuring more precisely

計測要素技術の開発
Development of
measurement element
technology

先端計測分析システムの開発
Development of advanced
measurement and analysis
system

新規開発
New development
ニーズを反映した改良
Improvements that reflect your needs

分析技術に仕上げる
Completion of analytical technology
分析評価
適用・応用
Analysis evaluation and application

公開利用・適用
Common use and research
support

大学や企業の研究開発、材料開発等を支援
Supporting academic and commercial R&D
and material development

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ANCF Equipment and Analysis Targets

Positron probe microanalyzer (PPMA)	nano vacancy	P.6
X-ray absorption fine structure spectroscopy with a superconducting fluorescence detector (SC-XAFS)	nano local structure, electronic state	P.10
Visible/near-infrared transient absorption spectrometer (VITA)	dynamics of photo-reaction	P.14
Real surface probe microscope (RSPM)	nano surface structure	P.18
Solid-state nuclear magnetic resonance spectrometer (SSNMR)	local structure of solid-state materials	P.22
Photoelectron spectroscopy with a pulse EUV source (EUPS)	nano surface physical properties	P.26
Superconducting scanning electron microscope (SC-SEM)	state analysis for light elements	P.30

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PPMA is a scanning positron microscope for mapping positron lifetimes in samples. The PPMA uses an intense positron microbeam (beam diameter $< 100 \mu\text{m}$) generated by an electron linear accelerator.

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
Positron Probe MicroAnalyzer (PPMA)

An intense positron beam generated by an electron linear accelerator is formed into a series of short pulses which are focused to micrometer size, and are then directed onto a sample. The structures of atomic and nanoscale vacancy-type defects and voids are evaluated by employing positron annihilation lifetime spectroscopy (PALS).

Time resolution	200 ps - 300 ps
Implantation energy	1 keV - 30 keV
Beam diameter	0.01 mm - 10 mm
Pulse duration	$\cong 30 \text{ ns}$

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


Schematic of positron lifetime spectroscopy (crystalline materials)

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

• Evaluation of relationship between hydrogen embrittlement and defects in high-strength steels³

The relationship between hydrogen embrittlement and defects in high-strength steels was investigated using the positron probe microanalyzer. Since the positron lifetime became longer in the vicinity of fracture surfaces for the tensile-stressed and hydrogen-charged specimen, it was found that vacancy-type defects and defect clusters are formed in the specimen.

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

Stretched polyethylene specimen

• Evaluation of structural change of stretched polymers⁴

The positron lifetimes of stretched low- and high-density polyethylenes were measured at different positions. It was found that the lifetime of the ortho-positronium (a bound state of an electron and a positron with parallel spins), which annihilates in amorphous regions of the polyethylene, was shortened in the stretched part. This indicates that the state of free volumes in the amorphous region changes due to stretching.

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
Variation of o-Ps lifetime upon stretching (LDPE)

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• Structural analysis of a-Si:H/c-Si heterojunction solar cell⁵

The size of micro voids in a-Si:H was measured using the positron probe microanalyzer, and the relationship between the void size and an optical parameter was investigated. It was found that the void size as determined by the positron lifetime changes depending on the film formation temperature, and that the size is correlated with the optical parameter.

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Relationship between void size in a-Si:H and optical parameter

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- フラーレン粉末体の空隙評価
- 構造材料の塑性変形、疲労破壊過程における空孔形成メカニズムの解明
- 逆浸透膜の空隙評価
- 酸化物半導体薄膜の欠陥解析
- Evaluation of pores in fullerene powders
- Elucidation of the mechanism of formation of vacancy-type defects and voids in the plastic deformation and fatigue fracture processes of structural materials
- Evaluation of free volumes in reverse osmosis membranes
- Defect analysis of oxide semiconductor films

参考文献 / References

- ¹ 産総研 TODAY, 14-9, (2014).
- ² 大島永康ら, 陽電子科学会会報, 陽電子マイクロプローブ, **1**, 27-39 (2013).
- ³ T. Doshida *et al.*, ISIJ International, Enhanced Lattice Defect Formation Associated with Hydrogen and Hydrogen Embrittlement under Elastic Stress of a Tempered Martensitic Steel, **52 (2)**, 198-207 (2012).
- ⁴ T. Oka *et al.*, Appl. Phys. Lett., Free volume change of elongated polyethylene films studied using a positron probe microanalyzer, **101**, 203108 (2012).
- ⁵ N. Matsuki *et al.*, ECS Transactions, A Novel Optical Characterization of a-Si:H/c-Si Interface Microstructures Based on Data of Positron Annihilation Spectroscopy, **92(9)**, 21-24 (2019).

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SC-XAFS can measure the nanostructure of a light element dopant in a solid sample due to the high flux of synchrotron radiation as X-ray source and precise energy resolution of the superconducting detector.

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X-ray Absorption Fine Structure Spectroscopy with a Superconducting Fluorescence Detector (SC-XAFS)

The apparatus provides partial fluorescence yield X-ray absorption spectroscopy using a synchrotron light source and superconducting tunnel junction X-ray spectrometer, achieving an excellent X-ray energy resolution of 10–20 eV.¹

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
SC-XAFS apparatus mounted on KEK-PF BL-11A

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

X-ray Absorption Fine Structure Spectroscopy with a Superconducting Fluorescence Detector

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
Example of a fluorescent X-ray spectrum measured with a superconducting detector. For reference, a fluorescent X-ray spectrum measured by a silicon drift detector is also shown

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
Fluorescent X-ray spectrum of nitrogen doped SiC obtained with a superconducting detector. The characteristic X-ray of N-K is clearly separated from that of C-K

• Evaluation of SiC characteristics for practical use of high-voltage switch materials²

SiC is a compound semiconductor that is used as a high withstand voltage switch material. Nitrogen is doped into SiC by high-temperature ion implantation. It is necessary to reduce the resistance in order to suppress the power consumption. However, when nitrogen is doped at a high concentration, the activation rate decreases and the resistance value increases.

To clarify the location of nitrogen in the lattice after ion implantation and annealing, the fluorescence yield X-ray absorption spectrum at the nitrogen K-edge was measured.

Comparison with the simulation revealed that the doped nitrogen was replaced with carbon immediately after ion implantation, and that the structure was not changed by annealing.

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X-ray absorption spectrum of SiC at N-K edge

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
• Measurement of GaN to reduce the power consumption of luminescent materials

GaN is a compound semiconductor that is used as a light emitting material. In order to reduce power consumption, it is necessary to increase the activation ratio by doping with Mg etc., but the local structure of the doped Mg cannot be measured due to the interference of the Ga-L line adjacent to the Mg-K line.

The fluorescence X-ray spectrum of Mg-doped GaN was measured using the apparatus, and it was found that Ga-L and Mg-K rays could be clearly distinguished. The elimination of Ga-L line interference is expected to help elucidate the local structure of Mg and reduce the power consumption.

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VITA is a spectrometer to measure a transient absorption spectrum and luminescence lifetime from the visible to infrared wavelength region.

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Nano/picosecond transient absorption spectrometer

Detection of transient absorption spectra and their time profile of excited state, reaction intermediate, and charge carrier produced by pulse laser irradiation

Nano/picosecond luminescence lifetime spectrometer

Detection of luminescence spectra and lifetime by pulse laser irradiation

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
Visible/near-Infrared Transient Absorption Spectrometer

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
Schematic diagram of transient absorption spectroscopy

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
Schematic diagram of photoluminescence lifetime spectroscopy

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
• Deterioration test of Si solar cell¹

The amount of conduction carriers produced in Si solar cells was determined from the time dependence of light absorbance by the carriers. A decrease of the absorbance at the surface region of the cells after degradation was confirmed by using surface sensitive excitation light at shorter wavelengths. This result indicates that the degradation of Si solar cells occurs at their surface.

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Transient absorption measurement of Si solar cell

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
Excitation wavelength dependence of transient absorption of Si solar cell before and after degradation test

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• Revealing the excited-state dynamics of organic electroluminescent materials²

Highly efficient thermally activated delayed fluorescence (TADF) materials form specific species in the excited state. Based on this information, a new mechanism of TADF is proposed.

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
Transient absorption spectral map of highly efficient TADF molecule (left) and new TADF mechanism (right)

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• Photoluminescence dynamics of room-temperature phosphorescent materials³

The photophysical properties of room-temperature phosphorescent (RTP) materials were investigated using a fluorescence lifetime spectrometer. For three kinds of benzil derivatives, mono benzil derivative 1 was revealed to exhibit the most efficient RTP.

室温燐光材料の化学構造式 (左) と室温燐光の時間プロフィール (右)
Room-temperature phosphorescent materials (left) and their luminescence time profiles (right)

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Nanoparticles, perovskite materials, photocatalysts, metal-oxide semiconductors, optical light fiber materials, nanocarbon materials, photosensitive materials, circular polarization materials, organometallic materials, cells, proteins, fluorescence probe materials

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Scanning probe microscope system with calibration system and sample preparation facility, such as reference material, laser interferometer and etc. The system includes functions, such as KPFM, elastic modulus, spectroscopic images in various environment (in air, in vacuum, and in liquid).

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
RSPM 1

High resolution nearfield scanning optical microscope. The optical resolution is 10nm in FTIR mode.

RSPM 2

Scanning probe microscope (SPM) system for nano material characterization (KPFM, SNDM, elastic modulus, and spectroscopic images) in various environment, such as in air, in liquid, and in vacuum.

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
Laser interferometer system for determining cantilever spring constant

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Reference material for calibrating carrier concentration in silicon

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

Real-time observation of developing process of photoresist

• Developing process of photoresist

The developing process of photoresist on Si substrate was observed in a liquid. The origin of defects in the electron beam lithography process was studied with a high-speed scanning system. Real-time AFM images were obtained every ten seconds in the developer solution. The photoresist, which is designed as a 50-nm line, was clearly observed as shown in the figures.

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Measurement of removal force of nano-particles with scanning probe microscope

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
• Removal force of nano-particles on Si substrate

The removal force of nano-particles on Si substrate was measured with a scanning probe microscope. Forces required to remove particles and nano-patterns were compared. The results will be useful for developing cleaning equipment for semiconductor factories.

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(Left) High resolution AFM images of CNT and DNA

(Right) Local electric-property measurement of semiconductor devices

- Surface morphology measurement of nano-materials such as CNT and DNA
- DNA or large molecules can be observed in liquid
- Surface potential and carrier concentration in silicon can be observed with KPFM, SCM and SSRM
- Nano-device measurement in contact mode (current, thermal property) is possible
- High-voltage surface charge can be measured with KPFM
- Nano-devices can be measured while in operation with a device analyzer

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SSNMR is used to measure the local structure and motion of atoms and molecules in solid materials by using atomic nuclear spin as a probe.

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Equipment 1


Solid-state NMR (600 MHz) for ^1H - ^{109}Ag , Wide bore SCM type, 1D and 2D experiments

Equipment 2

Solid-state NMR (200 MHz) for ^1H - ^{15}N , Wide bore SCM type, 1D and 2D experiments

Equipment 3


Solid-state NMR (20 MHz) for ^1H only, Permanent magnet type, relaxation time measurement

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
Solid-State Nuclear Magnetic Resonance Spectrometer

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
Magic Angle Spinning takes an average of phenomena such as dipole-dipole interactions, so a narrower spectrum will appear.

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In high-resolution solid-state NMR experiments, a special small sample tube (ZrO₂) is used, and the maximum sample rotation speed depends on the diameter of the tube selected.

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
A table of nuclides in "Solid-State NMR Spectral Database". It contains multi-nuclear NMR spectra of solid samples measured in AIST.

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Molecular structure of organic di-silane compound and part of the ^{13}C SSNMR spectrum

• Structure of light emitting materials (^{13}C)

Organic di-silane molecules are synthesized, and show a light emitting property and take various molecular conformations in the solution state. We have investigated molecular conformations in the solid state using solid-state NMR measurements. ^{13}C CP/MAS NMR was performed at 150.97 MHz using a 4 mm sample tube and 8 kHz rotation speed.


The figure shows the molecular structure of an organic di-silane compound and the temperature dependence and part of the ^{13}C SSNMR spectrum. This shows that no conformation change occurs in solid-state samples.

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^{31}P MAS NMR spectrum of BN with TMPO after ball-milling process

• Acid-base property of BN catalyst (^{31}P)

Hexagonal boron nitride (BN) exhibits nitro-aldol reaction and Knoevenagel condensation reaction activity after the ball-milling process. Absorbing basic probe molecules, trimethylphosphineoxide (TMPO), ^{31}P MAS NMR was performed at 161.98 MHz using a 4 mm sample tube and 8 kHz rotation speed.


The figure shows the ^{31}P MAS NMR spectrum: 63 ppm and 43 ppm signals were observed. The 43 ppm signal is assigned to physisorbed TMPO and the 63 ppm signal is assigned to adsorption on acidic sites, based on chemical shift values.

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^{13}C CP/MAS NMR spectra of SiPs, SiPs-initiator and SiPs-initiator-HMDS

・ Organic-inorganic hybrid materials (^{13}C)

The development of grafted silica particles of electrically conductive polymers was investigated. The figure shows the solid-state ^{13}C CP/MAS NMR spectrum of silica particles in preparing grafted silica particle components, SiPs, SiPs-initiator and SiPs-initiator-HMDS. Capping of -OH was achieved by replacing -HMDS residue.

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- 金属水素化物の構造解析評価、水素貯蔵材料における水素のサイトと拡散挙動
- 無機有機ハイブリッド材料における界面の構造、可逆的な溶解固化反応の構造変化を測定
- 電池材料の伝導メカニズムを評価、イオン伝導材料におけるイオン拡散とそのメカニズム
- 食肉の固体脂含量を評価
- Structural analysis of metal hydrate, hydrogen sites and the diffusion process in hydrogen storage materials
- Interfacial structure and change of structure in the reversible melt-solidification reaction in inorganic and organic hybrid materials
- Mechanisms of conduction of batteries, ionic diffusion, and ion conducting materials
- Quantification of solid fat included in meats

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EUPS is a photoelectron spectroscopy system, in which an EUV pulse (255.17 eV, 3-ns duration) emitted from a laser-produced plasma excites a sample and the energy spectrum of electrons is analyzed with a time-of-flight (TOF) analyzer. The principle of EUPS was invented and many technologies for implementing the EUPS were developed at AIST.¹ (This equipment is publicly open only under AIST-self-sponsored project.)

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

Photoelectron Spectroscopy with a Pulsed EUV Source (EUPS)

A high-power short-pulse laser is focused on a BN round bar target to generate a high-temperature plasma. The emission line at 255.17 eV from boron ions in the plasma is focused on a sample with an elliptical mirror. The electrons emitted from the sample travel in a flight tube, and the time taken to arrive at an electron detector is converted into kinetic energy to obtain the energy spectrum.

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Sample holder

4 samples (10 mm × 3 mm) can be loaded. Powder samples can be measured. Maximum size of sample: 20 mm × 20 mm. Temperature dependence from RT to 400 °C. Pre-treatment: sputtering, heating up to 1,000 °C in gas atmosphere

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

Photoelectron Spectroscopy with a Pulsed EUV Source

原理、特徴 / Principle and Features

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Schematic diagram of photoelectron spectroscopy

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
Relation of electron mean free path and electron energy

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
Photoelectron spectroscopy, which obtains the energy spectrum of photo-emitted electrons, provides information on the electronic states of the surface.


- In EUPS, the photon energy of the excitation light is 250 eV and the kinetic energy of the photoelectrons is in the energy region where the electron escape depth is the shallowest at 0.5 nm, so it has ultra-high sensitivity for detecting impurities of one atomic layer or less.
- Utilizing the pulse nature of the light source, the electron energy is analyzed by the time-of-flight (TOF) method. Since electrons in the entire energy region are measured at the same time, a wide spectrum can be obtained with a small number of irradiations.
- Thanks to this feature, the charging of insulators is minimized, and so the work function of insulators can be measured.
- Since the brightness of the light source is extremely high and the instantaneous photon density on the surface is very high, the resistance of an ultra-thin film can be evaluated. Band bending and charge trap centers having a lifetime of milliseconds or less can also be evaluated from the dependence of the electron spectrum on the excitation intensity.

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
In TOF, a spectrum in a wide energy region is obtained with only one-shot EUV irradiation.

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Detection of changes of the properties of TiO₂ nanocrystals by ozone cleaning of the Si substrate

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4 samples: Porous and compact TiO₂ with and without ozone cleaning of the Si substrate

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
The shift of the cut-off position, as seen in the left figure, gives:


a) the resistivity of a sample from the slope, and b) the density of charge trap centers from the rapid saturation at low intensity. The right figure shows that, in nano-porous TiO₂, the number of charge trap centers was reduced by ozone cleaning.

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The spectral shape near the cut-off reflects the energy relaxation time of excited electrons. In porous TiO₂, ozone cleaning increased the lifetime, indicating that the density of defects was reduced.

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

EUPS has high surface sensitivity for detecting sub-mono atomic-layer contamination owing to a small photoelectron escape depth of about 0.5 nm. Detection by EUPS of a very small amount of contaminants on a protective oxide of a plasma display panel, which was difficult to measure with ordinary XPS, led to identification of the contamination source, and the manufacturing process was improved.

Sub-monolayer organic molecule on MgO

Removal of contamination on MgO by heating

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¹ Development of EUPS for analyzing electronic states of topmost atomic layer

—Materialization of laser-produced plasma source application and EUPS observed fascinating surface—

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Scanning electron microscope with a superconducting tunnel junction X-ray detector (SC-SEM) is a SEM-EDX analyzer utilizing the superconducting X-ray detector which is realized a high throughput (10^{-3} Sr) of semiconductor X-ray detectors and a high energy resolution (<10 eV) of wave-dispersive X-ray detectors (WDS). The SC-SEM is suitable to obtain spatial information on the elemental composition of a sample nondestructively.

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
Scanning Electron Microscope with a Superconducting tunnel junction X-ray detector(SC-SEM)

- X-ray detector: Superconducting tunnel junction X-ray detector + Polycapillary X-ray lens
- SEM: JSM-IT800SHLs (JEOL)

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
Scanning Electron Microscope with a Superconducting tunnel junction X-ray detector

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
Energy spectra in the vicinity of the N-K α peak acquired using the STJ detector (red solid), SDD (green dots), and WDS (blue dashes).

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
Energy spectra for PPS acquired using the STJ detector (red line) and Si(Li)-SSD (blue line).

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• Elemental analysis of boron nitride (BN)¹

X-ray spectrum acquired by using the SC-SEM for a BN powder is displayed in the figure. Separate B- $K\alpha$ (188 eV), C- $K\alpha$ (277 eV), N- $K\alpha$ (393 eV), and O- $K\alpha$ (525 eV) peaks can be clearly observed due to the high energy resolution. The C- $K\alpha$ and O- $K\alpha$ peaks originated from carbon compound materials and an oxidized layer on the surface of the BN, respectively.

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
• Elemental analysis of a lithium fluoride (LiF)³

X-ray spectrum acquired by using the SC-SEM for a LiF crystal is displayed in the figure. Li- $K\alpha$ (54 eV) peaks can be detected due to the high sensitivity for low energy. The energy resolution of Li- $K\alpha$ was approximately 8 eV.

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
● Elemental analysis and chemical state analysis of precipitates in heat resistant steel^{2,4}

Light elements (B, C, N, O ...) have a big influence on performances of structural materials, such as high-resistance steels. To understand the effects of the light elements in the steels and to improve those properties, it is very important to evaluate amounts of each light element and their spatial distribution. Fig. (a) shows a SEM image of a 9% Cr steel. Black and white inclusions of sub-micrometer size can be seen. To evaluate the elemental composition of the black inclusion, X-ray spectra of point 1 and 2 in the SEM image were measured by using the SC-SEM (Fig. (b)). In the spectrum of point 1, Cr-LI (500 eV) and O-K α (525 eV) peaks can be clearly separated. In contrast, in the spectrum of point 2, the N-K α peak was observed in addition to peaks observed in the spectrum of point 1, and the intensity of O-K α peak was about 3 times higher than that in the spectrum of point 1. As a result of comparison with the two spectra, it is supposed that the black inclusion is the steel containing a lot of oxygen and a little nitrogen.

(a) SEM 像と (b) 9% クロム鋼の X 線スペクトル
(a) SEM image and (b) X-ray spectra for the 9% Cr steel.

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参考文献 / References

- ¹ G. Fujii, et. al., X-ray spectrometry, 46, 325 (2017).
- ² G. Fujii, et. al., Journal of the Japan Society for Precision Engineering, 197 (2020).
- ³ M. Ukibe, et. al., Microscopy and Microanalysis, 24, 1036 (2018).
- ⁴ G. Fujii, et. al., Microscopy and Microanalysis, 23, 1088 (2017).
- ⁵ Ukibe, et. al., 75th Annual Meeting of The Japanese Society of Microscopy, abstract (2019).
- ⁶ T. Irimatsugawa, et. al., IEEE Transactions on Applied Superconductivity, 29 (2019).

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If your use or research results can be disclosed or if your data can be publicly shared, we will provide support under the program of Advanced Research Infrastructure for Materials and Nanotechnology (ARIM); if confidentiality is required, we will provide support as an independent operation of AIST.

産総研の自主運用
Self-sponsored projects of AIST

マテリアル先端リサーチインフラ事業
Advanced Research Infrastructure for
Materials and Nanotechnology
(ARIM) Program

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Common use

・利用者が機器を操作
User operates the equipment on own

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Collaborative study

・利用者と産総研が共同で実施
Jointly implemented by users and AIST staff

契約・約款に基づく支援 (有料)
Contract-based support (charged)

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STEP 1

会員登録とユーザー登録 Member registration and user registration

Web サイトから会員登録とユーザー登録の後、事前相談申請をお願い致します。(登録等は無料です)
After registering as a member and user on the website, please apply for prior consultation (registration is free).

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We will introduce the person in charge of the equipment; please ask us about the details and form of support.

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If it is judged that support is possible, please apply.

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After use, you will be billed (collection of usage fee) according to the usage time.

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For use under the ARIM, you will be required to submit a usage report, register data, and complete the questionnaire form after the support ends.

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We provide facility training, lab tours and seminars. Please join us.

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Please see the latest information on the **Web** !

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National Institute of Advanced Industrial
Science and Technology
AIST Nanocharacterization Facility (ANCF)
Secretariat

AIST Tsukuba Central 2, 1-1-1 Umezono,
Tsukuba, Ibaraki 305-8568 Japan
E-mail: ancf-contact-ml@aist.go.jp